

# Sip/Puff Switch

**No Batteries or Bluetooth Pairing**

## Getting Started:

1. Attach your interface of choice – head mount or flexible shaft to one end of the air tube.
2. If using the headset, attach the opposite end of the air tube to the in-line filter. Then attach the in-line filter to the Switch Box. **NOT USING AT LEAST ONE FILTER WILL VOID THE WARRANTY**
3. If using the flexible shaft attach the air tube directly to the Switch Box. *There is a filter in the mouthpiece.*
4. Plug a stereo 3.5mm cable into the center Sip/Puff Switch connector if your device (being accessed) has a stereo, or two switch interface using one cable. If your device has two single switch inputs, then plug a 3.5mm cable into both Sip/Puff Switch connectors.
5. Plug the other end of the cable(s) into the device being accessed.
6. Use the Sip and Puff test buttons on the Switch Box to verify connectivity to your device, without sipping or puffing on the air tube.

The Sip/Puff Switch Box contains two pressure switches that are designed to close two sets of contacts when sipping or puffing on the air tube. The Sip/Puff Switch does not use power and therefore does not have batteries, an external adapter, or use power from the device being accessed. It is completely passive.

The rarely used slide switch will allow you to change the sense of sip and puff on the stereo cable to your device.

The electrical switches are designed for switching no more than 30-volts and 40-milliamps.

## Warranty:

A new Sip/Puff Switch is warranted against defects in material and workmanship for one (1) year from the date of shipment from Origin Instruments Corporation (OIC) or one of our dealers. See our full Sip/Puff Switch User Guide for additional warranty terms.

## Application Disclaimer:

The Sip/Puff system shall not be used in an application where its failure would lead to injury or property loss.



## Disposal:

At end of life these products should be taken to a recycling center. There are a variety of plastic materials, metals and a printed circuit board. All materials conform with EU RoHS and REACH requirements.

The Sip/Puff Switch, Headset and flexible shaft should not be discarded in regular consumer waste. In some countries you may return the product to your retailer.



May 26, 2021

Origin Instruments Corporation  
854 Greenview Dr  
Grand Prairie, TX 75050, USA  
<https://www.ori.com>  
972-606-8740